Product information

Information about other products is available at: www.demeditec.com

Renin ELISA

DE5125

96

Demeditec Diagnostics GmbH
Lise-Meitner-Strasse 2
24145 Kiel – Germany
www.demeditec.com

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

| Introduced modifications / Durchgeführte Änderungen / Modifiche introdotte / Modificaciones introducidas / Modifications apportées / Modificações introduzidas | |
|--|---|
| <p>The following changes have been made in comparison to the previous version: Im Vergleich zur Vorgängerversion wurden folgende Änderungen vorgenommen: Rispetto alla versione precedente, sono state apportate le seguenti modifiche: Se han introducido los siguientes cambios en comparación con la versión anterior: Les modifications suivantes ont été apportées par rapport à la version précédente : Foram efetuadas as seguintes alterações em comparação com a versão anterior:</p> | |
| <p>Detailed editorial revision. Changed wording in several chapters. Ausführliche redaktionelle Überarbeitung. Geänderter Wortlaut in mehreren Kapiteln. Revisione editoriale dettagliata. Modificato il testo in diversi capitoli. Revisión editorial detallada. Se ha cambiado la redacción de algunos capítulos. Révision éditoriale détaillée. Modification de la formulation dans plusieurs chapitres. Revisão editorial detalhada. Texto alterado em vários capítulos.</p> | |
| 4.2 Materials required but not provided: | Microplate shaker, 400 rpm (old: 300 - 700 rpm); |
| 5.2 Samples Storage: | Addition of the following: "To inactivate activated Prorenin, please incubate sample over night at 37 °C and repeat measurement." |
| 5.3 Sample Preparation: | Important note added for the cause of values above the highest standard. |
| 6.2 Test Procedure: | In step 4 and 7 shaker speed changed to 400 rpm (old: 300-700 rpm). |
| 6.3.1 Example of Typical Standard Curve: | updated |
| 7 REFERENCE VALUES: | For the 2.5 th percentile and for the minimum value of the range, the value is updated to <2.5 pg/mL, according to section 9.3. |
| 9.3 Sensitivity: | updated, LoB (2.499 pg/mL) is now given as lowest value for the measuring range, instead of previous listed analytical sensitivity (0.80 pg/mL). |
| 10.2 Drug Interferences: | updated |
| 9.6 Method Comparison: | added |
| 12 LITERATURE: | addition of ref. # 14. |

Table of Contents / Inhaltsverzeichnis / Tabella dei Contenuti / Tabla de Contenidos / Table des matières

| | | | | | |
|----|---|----|----|--|----|
| 1 | INTENDED USE | 4 | 8 | CONTROLLO DI QUALITÀ | 31 |
| 2 | PRINCIPLE OF THE TEST | 4 | 9 | CARATTERISTICHE PRESTAZIONALI | 32 |
| 3 | WARNINGS AND PRECAUTIONS..... | 5 | 10 | LIMITAZIONI DELLA PROCEDURA | 32 |
| 4 | MATERIALS | 6 | 11 | ASPETTI LEGALI | 33 |
| 5 | SAMPLE COLLECTION, STORAGE AND PREPARATION | 7 | 1 | FINALIDAD PREVISTA | 34 |
| 6 | ASSAY PROCEDURE | 8 | 2 | PRINCIPIO DEL TEST | 34 |
| 7 | REFERENCE VALUES | 10 | 3 | ADVERTENCIAS Y PRECAUCIONES..... | 34 |
| 8 | QUALITY CONTROL..... | 11 | 4 | MATERIAL..... | 35 |
| 9 | PERFORMANCE CHARACTERISTICS | 11 | 5 | TOMA DE LA MUESTRA, ALMACENAMIENTO Y PREPARACIÓN..... | 37 |
| 10 | LIMITATIONS OF THE PROCEDURE | 13 | 6 | PROCEDIMIENTO DEL ENSAYO..... | 38 |
| 11 | LEGAL ASPECTS | 14 | 7 | VALORES DE REFERENCIA..... | 40 |
| 1 | ZWECKBESTIMMUNG | 15 | 8 | CONTROL DE CALIDAD..... | 40 |
| 2 | TESTPRINZIP | 15 | 9 | CARACTERÍSTICAS DE FUNCIONAMIENTO | 41 |
| 3 | WARNUNGEN UND VORSICHTSMAßNAHMEN..... | 15 | 10 | LIMITACIONES DEL PROCEDIMIENTO | 41 |
| 4 | MATERIALIEN..... | 16 | 11 | CUESTIONES LEGALES | 42 |
| 5 | ENTNAHME, LAGERUNG UND VORBEREITUNG DER PROBEN | 19 | 1 | DESTINATION DU DISPOSITIF..... | 43 |
| 6 | TESTDURCHFÜHRUNG..... | 20 | 2 | PRINCIPE DU TEST | 43 |
| 7 | REFERENZWERTE | 22 | 3 | AVERTISSEMENTS ET PRECAUTIONS..... | 43 |
| 8 | QUALITÄTSKONTROLLE | 23 | 4 | MATÉRIAUX..... | 44 |
| 9 | LEISTUNGSMERKMALE | 23 | 5 | PRÉLÈVEMENT, STOCKAGE ET PRÉPARATION DES ÉCHANTILLONS | 46 |
| 10 | GRENZEN DES VERFAHRENS | 23 | 6 | PROCÉDURE DE DOSAGE | 47 |
| 11 | RECHTLICHE GRUNDLAGEN..... | 24 | 7 | VALEURS DE RÉFÉRENCE..... | 49 |
| 1 | DESTINAZIONE D'USO | 25 | 8 | CONTRÔLE DE QUALITÉ | 49 |
| 2 | PRINCIPIO DEL TEST | 25 | 9 | CARACTERISTIQUES EN MATIERE DE PERFORMANCES | 50 |
| 3 | AVVERTENZE E PRECAUZIONI..... | 25 | 10 | LIMITES DE LA PROCÉDURE..... | 50 |
| 4 | MATERIALI..... | 26 | 11 | ASPECTS JURIDIQUES | 51 |
| 5 | PRELIEVO, CONSERVAZIONE E PREPARAZIONE DEI CAMPIONI..... | 28 | 12 | LITERATURE | 52 |
| 6 | PROCEDURA DEL DOSAGGIO | 29 | | SYMBOLS USED WITH DEMEDITEC ASSAYS... | 52 |
| 7 | VALORI DI RIFERIMENTO | 31 | | | |

1 INTENDED USE

The **DEMEDITEC Renin ELISA** is a manual enzyme immunoassay for the **quantitative** measurement of active Renin in human serum or EDTA plasma.

For *in vitro* diagnostic use. For laboratory professional use.

1.1 Scientific Validity Report

Renin is an enzyme (Mw of 37 kDa) that belongs to the aspartic acid protease family. Renin is a member of Renin-Angiotensin-Aldosterone System (RAAS) that controls blood pressure, renal blood flux, glomerular filtration, and sodium/potassium homeostasis. Renin is produced constitutively as prorenin, an inactive precursor with 386 amino acids, in the juxtaglomerular cells of the kidney (1). In response to low intra-renal blood pressure, reduced sodium reabsorption, hypokalemia or activity of the sympathetic nervous system, active renin can be released either from a depot in the kidney or generated from prorenin by cleavage of 46 amino acids at the N-terminus of prorenin (2,3). Prorenin secretion into the blood is continuous, in contrast to the tightly controlled release of renin, and blood concentration of prorenin is approx. 100-fold higher than active renin (4,5). After release and activation, soluble renin mediates cleavage of the α 2-globulin angiotensinogen into the precursor peptide angiotensin I, which ultimately is processed by angiotensin converting enzyme (ACE) to the octapeptide angiotensin II. All actions of angiotensin II are mediated by the G protein-coupled angiotensin type 1 (AT1) and angiotensin type 2 (AT2) receptors (6). Direct physiological effects of Angiotensin II include vasoconstriction, increase of tubular reabsorption of sodium and chloride, water retention, and release of the hormones aldosterone from adrenal cortex, antidiuretic hormone (ADH, Vasopressin) from posterior pituitary, and adrenocorticotrophic hormone (ACTH, Corticotropin) from anterior pituitary. Release of these hormones further supports sodium retention and secretion of potassium/H⁺ in the kidney, and increases thirst sensation and the desire for salt through the subfornical organ of the brain (7,8). In a negative feedback loop, renin secretion is reduced by high concentration of angiotensin II (9), and release of aldosterone is lowered by potassium depletion (10). Beside the action of soluble renin, binding of renin and prorenin to the membrane-bound renin receptor ATP6AP2 in brain, heart, placenta, liver, kidney and pancreas enhances efficiency of angiotensinogen cleavage and induces angiotensin-independent intra-cellular effects by activating mitogen activated kinases ERK1 and ERK2 (11). Besides the effects of the renin-angiotensin system (RAS) in the field of cardiovascular physiology as the main player in blood pressure homeostasis, other effects have since been described, and include proliferation, fibrosis, and inflammation (14).

Plasma renin is a good index for the activity of the RAAS. In case of dysfunction of the RAAS, the Renin assay will allow clinical implications for diagnosis, treatment, and follow up. Active renin should be measured in:

- Diagnosis of hypertension (high blood pressure: if diastolic blood pressure is > 90 mm Hg and systolic blood pressure is > 140 mm Hg; guideline of the European Society of Cardiology and the European Society of Hypertension)
- Differential diagnosis of hyperaldosteronism (primary hyperaldosteronism, secondary hyperaldosteronism with or without hypertension, pseudo-hyperaldosteronism)
- Diagnosis of isolated deficit in mineral corticoids
- Differential diagnosis of hypokalemia (secondary hyperaldosteronism or primary hypermineralcorticoism)
- Detection of Renin producing tumors in the kidney
- Monitoring of glucocorticoid therapy
- Diagnosis of insufficient response to antihypertensive treatment

2 PRINCIPLE OF THE TEST

The DEMEDITEC Renin ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **sandwich principle**. The microtiter wells are coated with a monoclonal antibody (mouse) directed towards a unique antigenic site of the human active Renin molecule. During the first incubation (together with assay buffer), Renin in the added sample binds to the immobilized antibody. After incubation, unbound components are washed off. The afterwards added enzyme conjugate, which contains monoclonal anti-Renin antibody conjugated to horseradish peroxidase, binds to the Renin forming a sandwich complex.

After a washing step to remove all unbound enzyme conjugate, the solid phase is incubated with the substrate solution. The colorimetric reaction is stopped by addition of stop solution, and optical density (OD) of the resulting yellow product is measured. The intensity of color is proportional to the concentration of the analyte in the sample. A standard curve is constructed by plotting OD values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

3 WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* diagnostic use only. For laboratory professional use only.
- Before starting the assay, read the instructions for use completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
- Do not mix or use components from kits with different lot numbers. It is advised not to interchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Do not use reagents beyond expiry date as shown on the kit labels.
- Do not reuse microtiter wells.
- Reagents of other manufacturers must not be used together with the reagents of this test kit.
- All reagents in this kit are clear liquids, substrate solution is clear and colorless. Changes in its appearance may affect the performance of the test. In that case, contact DEMEDITEC.
- Microbial contamination of reagents or samples may give false results.
- Allow the reagents to reach room temperature (20 °C to 26 °C) before starting the test. Temperature will affect the optical density readings of the assay.
- All indicated volumes must be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution coloured. Do not pour reagents back into original vials as reagent contamination may occur.

General precautions

- Follow laboratory quality assurance and laboratory safety guidelines.
- Never pipet by mouth and avoid contact of reagents and samples with skin and mucous membranes.
- Do not smoke, eat, drink, or apply cosmetics in areas where samples or kit reagents are handled.
- Wear lab coats and disposable latex gloves when handling samples and reagents and where necessary safety glasses.

Biohazard information

- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. However, no known test method can offer total assurance that no infectious agent is present.
- The device contains material of animal origin, which is certified apparently free of infectious or contagious diseases and injurious parasites.
- Bovine components originate from countries where BSE (Bovine spongiform encephalopathy) has not been reported.
- All materials and samples of human or animal origin must be handled as if capable of transmitting infectious diseases.
- Handling must be done in accordance with the procedures defined by appropriate national biohazard and safety guideline or regulation. Waste must be discarded according to local rules and regulations.

Information to chemical hazards and hazard classification

- Some reagents contain preservatives in non-declarable concentrations. Nevertheless, in case of contact with eyes or skin, flush immediately with water.
- Substrate Solution contains an ingredient in non-declarable concentrations which causes serious eye irritation. In case of possible contact with eyes, rinse immediately carefully and thoroughly with eye wash or water. After contact with skin, wash with plenty of water. Take-off contaminated clothing and wash it before reuse.

- Avoid contact with Stop Solution containing < 5 % H₂SO₄. It may cause skin irritation and burns.
- Chemicals and prepared or used reagents must be treated as hazardous waste according to the national safety guideline or regulation.
- This product does not contain substances which have carcinogenic, mutagenic or toxic for reproduction (CMR) properties.

4 MATERIALS

4.1 Materials provided with the kit

1. **SORB MT Microtiterwells**, 12 x 8 wells (break apart), ready to use; Microtiter plate coated with anti-Renin antibody (monoclonal).
2. **CAL 0 – 5 LYO Standard (Standard 0-5)**, 6 x 1 mL, lyophilized; See “Reagent Preparation”, Standards*, concentrations: 0 – 4 – 16 – 32 – 64 – 128 pg/mL; Conversion: 1 pg/mL = 1.44 µIU/mL; *Calibrated against the following reference material: WHO 1st International Standard for Renin 68/356.*
3. **CONTROL low & high LYO Control Low & Control High**, 2 x 1 mL, lyophilized; See „Reagent Preparation“. Controls*; *For control values and ranges please refer to vial label or Certificate of Analysis.*
4. **BUF Assay Buffer**, 1 x 20 mL, ready to use; Assay Buffer*.
5. **ENZ CONJ Enzyme Conjugate**, 1 x 14 mL, ready to use; Enzyme Conjugate*, monoclonal anti-human Renin antibody conjugated to horseradish peroxidase.
6. **SUB TMB Substrate Solution**, 1 x 14 mL, ready to use; Substrate Solution, Contains 3,3',5,5'-tetramethylbenzidine (TMB). *Keep away from direct sun light.*
7. **STOP SOLN Stop Solution**, 1 x 14 mL, ready to use; Stop Solution, contains < 5 % H₂SO₄, *avoid contact with the stop solution. It may cause skin irritations and burns.*
8. **WASH SOLN 40x Wash Solution**, 1 x 30 mL, See “Reagent Preparation”; Wash Solution, 40X concentrate.
9. 1x **Instructions for Use**; 1x **Certificate of Analysis (CoA)**

* Contains non-mercury preservative.

4.2 Materials required but not provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm)
- Calibrated variable precision micropipettes
- Incubator for 37 °C
- Manual or automatic equipment for rinsing microtiter plate wells
- Absorbent paper
- Distilled water
- Timer
- Graph paper or software for data reduction
- Microplate shaker, 400 rpm

4.3 Storage and Stability of the Kit

Unopened kits and reagents as well as **opened reagents** must be stored at 2 °C to 8 °C.

The microplate must always be stored in the resealable aluminum pouch containing a desiccant. Do not open the pouch until it has reached room temperature. The microtiter plate consists of 12 individual strips. Each strip can be divided into 8 individual wells. Unused wells must be immediately returned to the aluminum pouch with the desiccant and stored again tightly resealed at 2 °C to 8 °C.

Once opened, reagent vials must be closed tightly again.

| | Storage Temperature | Stability |
|-------------------------------------|---------------------|--|
| Unopened kits and unopened reagents | 2 °C to 8 °C | Until the expiration date printed on the label. Do not use reagents beyond this date! |
| Opened kit | 2 °C to 8 °C | 8 weeks (For reconstituted reagents refer to “4.4 Reagent Preparation”.) |

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature (20 °C to 26 °C) prior to use.

Standards

Reconstitute the lyophilized contents of each standard vial with 1 mL distilled water and let stand for at least 10 minutes at room temperature. Mix several times before use.

| | | |
|---------------------------------|-------------------------|-----------------|
| Stability after reconstitution: | at 2 °C to 8 °C | 14 days |
| | at -20 °C (in aliquots) | up to 12 months |

Controls

Reconstitute the lyophilized content of each control vial with 1 mL distilled water and let stand for at least 10 minutes at room temperature. Mix several times before use.

| | | |
|---------------------------------|-------------------------|-----------------|
| Stability after reconstitution: | at 2 °C to 8 °C | 14 days |
| | at -20 °C (in aliquots) | up to 12 months |

Wash Solution

Add distilled water to the 40X concentrated Wash Solution. Dilute 30 mL of concentrated *Wash Solution* with 1170 mL distilled water to a final volume of 1200 mL.

| | | |
|---------------------------|-------------------|--------|
| Stability after dilution: | at 20 °C to 26 °C | 1 week |
|---------------------------|-------------------|--------|

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DEMEDITEC must be informed in writing, at the latest one week after receiving the kit. Damaged single components must not be used for a test run. They have to be stored until a final solution has been found. After this, they must be disposed of according to the official regulations.

5 SAMPLE COLLECTION, STORAGE AND PREPARATION

The following sample material can be used in this test:

Human serum or EDTA plasma

Samples containing sodium azide should not be used in the assay. In general, it should be avoided to use hemolytic, icteric, or lipemic samples. For further information, refer to chapter "*Interfering Substances*".

Conditions under which samples are collected must be carefully controlled, since a number of physiological factors can influence the renin secretion. These include:

- **Posture:** the patient must have been lying down for more than 1 hour or upright for more than 1 hour
- **Daily Renin oscillations:** sampling is to be done between 7 AM and 10 AM if possible.
- **Diet:** sodium content in the diet must be known and eventually verified by the measurement of natriuria over a period of 24 hours
- **Medication:** the level of active renin can be affected by antihypertensive medication (e.g. diuretics, ACE inhibitors, beta adrenergic blocking agents, vasodilators, renin inhibitors)
- **Pregnancy:** the level of inactive and active renin increases during pregnancy
- **Menstrual cycle:** the level of active renin increases on the second phase of the cycle (sampling is to be done if possible, during the first phase)
- **Age:** active renin level decreases with age

5.1 Sample Collection

Serum: Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma: Whole blood should be collected into centrifuge tubes containing anticoagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

Whole blood should not be frozen before centrifugation.

5.2 Samples Storage

Samples must be stored tightly capped prior to performing the assay. If stored frozen, freeze only once. Thawed samples must be inverted several times prior to testing.

| | | | |
|-----------|---------------------------|-----------------|--|
| Stability | at 2 °C to 8 °C | NOT stable | Store at room temperature and NOT at 2 °C to 8 °C prior to processing, since cryoactivation of prorenin may occur in the temperature range of 2 °C to 8 °C, giving false positive active renin values > 128 pg/mL (12,13, 14). To inactivate activated Prorenin, please incubate sample over night at 37 °C and repeat measurement. If samples cannot be tested <u>within 4 hours</u> of collection, store frozen. It is recommended to rapidly freeze, and thaw samples avoiding the temperature range of 2 °C to 8 °C. A dry ice/ethanol bath can be used for rapid freezing procedures. |
| | at room temperature | < 4 hours | |
| | at ≤ -20 °C (in aliquots) | up to 12 months | |

5.3 Sample Preparation

Samples can be assayed without additional preparation.

If in an initial assay, a sample is found to contain more analyte than the highest standard, the sample can be diluted with *Assay Buffer* and reassayed as described in the assay procedure.

Important note: values above highest standard (> 128 pg/mL) may be caused by:

1. very high concentrations of Renin (only in rare cases). In this case, please follow the instructions for dilution given below.
2. Cryoactivation of Prorenin to activated Prorenin due to inappropriate sample storage. In this case follow instruction given in chapter 5.2

If a dilution is required, the sample must be diluted at least 1:10 with *Assay Buffer*. For calculation of concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:10: 20 µL sample + 180 µL *Assay Buffer* (mix thoroughly)
- b) dilution 1:20: 10 µL sample + 190 µL *Assay Buffer* (mix thoroughly).

6 ASSAY PROCEDURE

6.1 Procedural Notes

- All reagents and samples must be allowed to come to room temperature (20 °C to 26 °C) before use.
- All reagents must be mixed without foaming.
- Do not interchange caps of reagent vials to avoid cross-contamination.
- Use new disposal plastic pipette tips for each standard, control, or sample in order to avoid carry-over.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Mix the contents of the microtiter plate wells thoroughly to ensure good test results.
- Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
- Once the test has been started, all steps must be completed without interruption and in the same sequence for each step.
- The enzymatic reaction is linearly proportional to time and temperature.

- Optical density is a function of the incubation time and temperature. Respect the incubations times and temperatures as given in chapter “Test Procedure”.
- Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- **Important note to wash procedure:**
Washing is critical. Improperly washed wells will give erroneous results. The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
- **Test performance using fully automated analysis devices:**
Automated test performance using fully automated, open-system analysis devices is possible. However, the combination must be validated by the user.

6.2 Test Procedure

Each run must include a standard curve. The controls serve as internal controls for the reliability of the test procedure. They must be assayed with each test run. The given test procedure describes manual processing.

Important note: The accuracy of this assay is markedly influenced by the correct incubation temperature, incubation time, and correct pipetting volumes.

1. Secure the desired number of microtiter wells in the frame holder.
2. Pipette **150 µL** of **Assay Buffer** in all appropriate wells.
3. Pipette **50 µL** of each **Standard, Control,** and **sample** with new disposable tips into appropriate wells.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
4. Incubate for **90 minutes** at room temperature on a plate shaker with 400 rpm.
5. Wash the wells as follows:
If the wash step is performed manually:
 Briskly shake out the contents of the wells.
 Rinse the wells **4 times** with **300 µL** diluted *Wash Solution* per well.
If an automated plate washer is used:
 Rinse the wells **4 times** with **400 µL** diluted *Wash Solution* per well. .
At the end of the washing step, always strike the wells sharply on absorbent paper to remove residual droplets!
6. Add **100 µL Enzyme Conjugate** into each well.
7. Incubate for **90 minutes** at room temperature on a plate shaker with 400 rpm.
8. Wash as described in step 5.
9. Pipette **100 µL** of **Substrate Solution** to each well.
10. Incubate for **15 minutes** at room temperature.
11. Stop the enzymatic reaction by adding **100 µL** of **Stop Solution** to each well.
12. Measure the optical density (OD) of the solution in each well at **450 nm (reading) and at 620 nm to 630 nm (background subtraction, recommended)** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

6.3 Calculation of Results

1. The concentration of the samples can be read directly from the standard curve.
2. For duplicate determinations, the mean of the two optical density (OD) values for each standard, control, and patient sample must be taken. If the two values deviate substantially from one another, DEMEDITEC recommends retesting the samples.
3. Samples with concentrations exceeding the highest standard can be further diluted with *Assay Buffer* and re-assayed as described in "Test Procedure", or must be reported as > 128 pg/mL. For the calculation of the concentrations, this dilution factor must be considered.
4. Automated method:
The results in the instructions for use have been calculated automatically using a four-parameter logistic (4PL) curve fit. (4PL Rodbard or 4PL Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. Manual method:
Using linear graph paper, construct a standard curve by plotting the (mean) OD obtained from each standard against its concentration with OD value on the vertical (Y) axis and concentration on the horizontal (X) axis.
Determine the corresponding sample concentration from the standard curve by using the (mean) OD value for each sample.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

| Standard | Optical Density (450 nm) |
|------------------------|--------------------------|
| Standard 0 (0 pg/mL) | 0.170 |
| Standard 1 (4 pg/mL) | 0.284 |
| Standard 2 (16 pg/mL) | 0.544 |
| Standard 3 (32 pg/mL) | 0.885 |
| Standard 4 (64 pg/mL) | 1.573 |
| Standard 5 (128 pg/mL) | 2.586 |

7 REFERENCE VALUES

It is strongly recommended that each laboratory determine its own reference values.

In a study conducted with apparently normal healthy adults, using the DEMEDITEC Renin ELISA the following values are observed in EDTA plasma:

| Healthy Adults | n | Mean (pg/mL) | Median (pg/mL) | 2.5 th - 97.5 th Percentile (pg/mL) | Range (min. - max.) (pg/mL) |
|------------------|----|--------------|----------------|---|-----------------------------|
| Supine position | 59 | 16.23 | 12.40 | < 2.50 – 53.83 | < 2.50 – 58.78 |
| Upright position | 59 | 19.73 | 16.18 | 2.79 – 61.83 | < 2.50 – 95.56 |

These values are also valid for serum.

In a study conducted with apparently normal healthy adults, using the DEMEDITEC Aldosterone ELISA (DE5298) and the DEMEDITEC Renin ELISA the following **Aldosterone-Renin Ratios** were determined in EDTA plasma:

Ratio Aldosterone Renin (pg/mL / pg/mL)

| | n | Mean | Median | 2.5 th - 97.5 th Percentile |
|----------------|----|------|--------|---|
| Healthy Adults | 89 | 8.68 | 5.30 | 0.52 – 37.83 |

These values are also valid for serum.

Values above or below the reference range should be considered as suspicious and require additional testing. The results alone should not be the only reason for any therapeutic consequences. The results must be correlated to other clinical observations and diagnostic tests.

8 QUALITY CONTROL

Good quality assurance in the laboratory requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results. Use controls at both normal and pathological levels. The controls and the corresponding results of the Quality Control Laboratory are stated in the Certificate of Analyses (CoA) added to the kit. The values and ranges stated on the CoA always refer to the current kit lot and must be used for direct comparison of the results. If available, it is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results. Apply appropriate statistical methods for analyzing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods. After checking the above-mentioned items without finding any error contact your distributor or DEMEDITEC directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Specificity of Antibodies (Cross-Reactivity)

The following substances were tested for cross-reactivity of the assay:

| Substance | Added concentration | Mean Cross-Reactivity (%) |
|---------------------|---------------------|---------------------------|
| Hepcidin | 140 ng/mL | < 0.001 |
| Pepsin | 50 µg/mL | < 0.001 |
| Human serum albumin | 1.2 mg/mL | < 0.001 |
| Human IgG | 1.0 mg/mL | < 0.001 |

Mean cross reactivity with Prorenin was 0.71% (mean value when prorenin was spiked in a concentration range from 256 pg/mL to 4096 pg/mL). However, the observed cross reactivity may only represent a contamination of the recombinant prorenin preparation with active renin due to auto-activation.

9.2 Sensitivity

| | |
|-------------------------------|-------------------------|
| Limit of Blank (LoB) | 2.499 pg/mL |
| Limit of Detection (LoD) | 4.308 pg/mL |
| Limit of Quantification (LoQ) | 6.021 pg/mL |
| Measuring range | 2.499 pg/mL – 128 pg/mL |
| Linear range | 2.850 pg/mL – 128 pg/mL |

9.3 Reproducibility

9.3.1 Within-run Precision

The within-run precision was determined with patient samples covering the complete measuring range within 20 days in 2 independent runs per day. CV was calculated as mean CV of 40 runs.

| Sample | n | Mean (pg/mL) | CV (%) |
|--------|----|--------------|--------|
| 1 | 40 | 106.03 | 1.5 |
| 2 | 40 | 71.01 | 1.5 |
| 3 | 40 | 48.52 | 2.0 |
| 4 | 40 | 30.13 | 2.4 |
| 5 | 40 | 20.27 | 1.9 |
| 6 | 40 | 10.62 | 4.8 |

9.3.2 Between-run Precision

The between-run precision was determined with patient samples covering the complete measuring range within 20 days in 2 independent runs per day and with 2 replicates per run (20 x 2 x 2). CV was calculated from 80 determinations.

| Sample | n | Mean (pg/mL) | CV (%) |
|--------|----|--------------|--------|
| 1 | 80 | 106.03 | 2.5 |
| 2 | 80 | 71.01 | 2.1 |
| 3 | 80 | 48.52 | 2.8 |
| 4 | 80 | 30.13 | 3.6 |
| 5 | 80 | 20.27 | 4.1 |
| 6 | 80 | 10.62 | 6.9 |

9.3.3 Between-lot Precision

The between-lot variation was determined by 6 measurements of different samples with 3 different kit lots.

| Sample | n | Mean (pg/mL) | CV (%) |
|--------|----|--------------|--------|
| 1 | 18 | 101.72 | 1.2 |
| 2 | 18 | 71.46 | 1.7 |
| 3 | 18 | 49.32 | 2.0 |
| 4 | 18 | 27.65 | 1.2 |
| 5 | 18 | 18.98 | 1.5 |
| 6 | 18 | 10.09 | 5.5 |

9.4 Recovery

Recovery was determined by adding increasing amounts of the analyte to different patient samples containing different amounts of endogenous analyte. The percentage recoveries were determined by comparing expected and measured values of the samples.

| | Sample 1 | Sample 2 | Sample 3 | Sample 4 |
|------------------------------|-------------|----------|----------|----------|
| Concentration (pg/mL) | 8.93 | 15.97 | 55.27 | 100.61 |
| Average Recovery (%) | 100.2 | 96.0 | 108.4 | 99.0 |
| Range of Recovery (%) | from | 95.5 | 86.8 | 106.2 |
| | to | 103.3 | 105.3 | 111.9 |

9.5 Linearity

Samples containing different amounts of analyte were serially diluted with *Assay Buffer*. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

| | Sample 1 | Sample 2 | Sample 3 | Sample 4 |
|------------------------------|-------------|----------|----------|----------|
| Concentration (pg/mL) | 45.16 | 53.20 | 67.35 | 126.02 |
| Average Recovery (%) | 101.7 | 102.8 | 110.0 | 98.5 |
| Range of Recovery (%) | from | 96.7 | 95.6 | 105.0 |
| | to | 108.6 | 114.6 | 114.6 |
| | | 100.8 | | |


9.6 Method Comparison

A comparison of DEMEDITEC Renin ELISA (DE5125) (y) and the reference method (x) (radioimmunoassay RENIN III GENERATION, Cisbio Bioassays) using clinical samples gave the following correlation:

$$n = 151$$

$$r = 0.979$$

$$y = 1.004x + 0.4047$$

10 LIMITATIONS OF THE PROCEDURE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the instructions for use and in compliance with the laboratory quality assurance guidelines. Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

Hemoglobin (up to 4 mg/mL), bilirubin (up to 0.5 mg/mL) and triglyceride (up to 30 mg/mL) have no influence on the assay results.

10.2 Drug Interferences

| Substance | Concentration Range of Spiked Substance | Mean Recovery (%) |
|-------------------------|---|-------------------|
| Aliskiren | 10 – 100 µg/mL | 139.9 |
| Captopril | 0.5 – 5 µg/mL | 102.4 |
| Furosemid | 6 – 60 µg/mL | 95.4 |
| Enalapril maleate | 0.05 – 0.5 µg/mL | 102.4 |
| Nicardipin-hydrochlorid | 5 – 50 µg/mL | 96.1 |
| Cathepsin D | 0.05 – 0.5 U/mL | 101.5 |
| Cathepsin D | 0.01 – 0.1 U/mL | 98.6 |

In addition, the level of active renin in plasma may be affected by antihypertensive medication (e.g. diuretics, ACE inhibitors, beta adrenergic blocking agents, or vasodilators)

10.3 High-Dose Hook Effect

"High-Dose Hook Effect" is not detected up to 8200 pg/mL of Renin.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the laboratory quality assurance guidelines and applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. If there is any doubt or concern regarding a result, please contact DEMEDITEC.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient. Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived. The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement. Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

11.4 Reporting of Serious Incident

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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
11.4 Notification des incidents graves

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12 LITERATURE

1. Imai T, Miyazaki H, Hirose S, et al. Cloning and sequence analysis of cDNA for human renin precursor. Proc. Natl. Acad. Sci. (1983) 80, 7405–7409.
2. Reudelhuber TL, Ramla D, Chiu L, et al. Proteolytic processing of human prorenin in renal and non-renal tissues. Kidney Int. (1994) 46, 1522–1524.
3. Neves FA, Duncan KG, Baxter JD. Cathepsin B is a prorenin processing enzyme. Hypertension (1996) 27, 514–517.
4. Müller DN, Luft FC. Direct renin inhibition with aliskiren in hypertension and target organ damage. Clin J Am Soc Nephrol. (2006) 1, 221-8.
5. Toffelmire EB, Slater K, Corvol P, et al. Response of plasma prorenin and active renin to chronic and acute alterations of renin secretion in normal humans. Studies using a direct immunoradiometric assay. J Clin Invest. (1989) 83, 679–687.
6. Carey RM, Padia SH. Angiotensin AT2 receptors: control of renal sodium excretion and blood pressure. Trends Endocrinol Metab. (2008) 19, 84-7.
7. Koepfen BM, Stanton BA. Renal Physiology (4th ed.). Philadelphia, PA. Mosby Physiology Monograph Series, 2007.
8. Navar LG, Inscho EW, Majid DSA, et al. Paracrine regulation of the renal microcirculation. Physiol. Rev. (1996) 76, 425–536.
9. Müller MW, Todorov V, Krämer BK, Kurtz A. Angiotensin II inhibits renin gene transcription via the protein kinase C pathway. Pflugers Arch. (2002) 444, 499-505.
10. Spät A, Hunyady L. Control of aldosterone secretion: a model for convergence in cellular signaling pathways. Physiol Rev. (2004) 84, 489-539.
11. Nguyen G., Delarue F., Burcklé C., et al. Pivotal role of the renin/prorenin receptor in angiotensin II production and cellular responses to renin. J Clin Invest. (2002) 109, 1417–1427.
12. Pitarresi TM., Rubattu S, Henrikson R, Sealey JE. Reversible cryoactivation of recombinant human prorenin. J.Biol.Chem. (1992) 267, 11753-9.
13. Nicar MJ. Specimen processing and renin activity in plasma. Clin. Chem. (1992) 38, 598.
14. Laghnam D. et al. Renin–Angiotensin–Aldosterone System and Immunomodulation: A State-of-the-Art Review. Cells (2021), 10, 1767.

SYMBOLS USED WITH DEMEDITEC ASSAYS

| Symbol | English | Deutsch | Française | Espanol | Italiano |
|---|------------------------------------|--|--|---|-------------------------------------|
|  | European Conformity | CE-Konformitätskennzeichnung | Conforme aux normes européennes | Conformidad europea | Conformità europea |
|  | Consult instructions for use | Gebrauchsanweisung beachten | Consulter les instructions d'utilisation | Consulte las Instrucciones | Consultare le istruzioni per l'uso |
|  | In vitro diagnostic device | In-vitro-Diagnostikum | utilisation Diagnostic in vitro | Diagnóstico in vitro | Per uso Diagnostica in vitro |
|  | For research use only | Nur Forschungszwecke für | Seulement dans le cadre de recherches | Sólo para uso en investigación | Solo a scopo di ricerca |
|  | Catalogue number | Katalog-Nr. | Référence | Número de catálogo | No. di catalogo |
|  | Lot. No. / Batch code | Chargen-Nr. | No. de lot | Número de lote | Lotto no |
|  | Contains sufficient for <n> tests/ | Ausreichend für "n" Ansätze | Contenu suffisant pour "n" tests | Contenido suficiente para <n> ensayos | Contenuto sufficiente per "n" saggi |
|  | Note warnings and precautions | Warnhinweise und Vorsichtsmaßnahmen beachten | Avertissements et mesures de précaution font attention | Tiene en cuenta advertencias y precauciones | Annoti avvisi e le precauzioni |
|  | Storage Temperature | Lagerungstemperatur | Température de conservation | Temperatura de conservación | Temperatura di conservazione |
|  | Expiration Date | Mindesthaltbarkeitsdatum | Date limite d'utilisation | Fecha de caducidad | Data di scadenza |
|  | Legal Manufacturer | Hersteller | Fabricant | Fabricante | Fabbricante |
| <i>Distributed by</i> | Distributed by | Vertrieb durch | Distribution par | Distribución por | Distribuzione da parte di |
|  | Version | Version | Version | Versión | Versione |
|  | Single-use | Einmalverwendung | À usage unique | Uso único | Uso una volta |