# Product information

Information about other products is available at: [www.demeditec.com](http://www.demeditec.com)

# Testosterone free ELISA

**DE2924**

**96**

Demeditec Diagnostics GmbH  
Lise-Meitner-Strasse 2  
24145 Kiel – Germany  
[www.demeditec.com](http://www.demeditec.com)

**CONTENTS / CONTENUTO / CONTENIDO**

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## 1. INTENDED PURPOSE

### **For *In Vitro* Diagnostic Use. For Laboratory Professional Use**

Free Testosterone ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of free testosterone in human serum or plasma from an adult population. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis and monitoring of disorders involving the male sex hormones (androgens).

## 2. CLINICAL SIGNIFICANCE

Testosterone is found in circulation predominantly linked to carrier proteins, the most common of which being sex-hormone binding globulin (SHBG). Testosterone plays a key role in the development of primary and secondary sexual characteristics in males and is involved in the production of female sexual hormones. Only 1 – 2% of testosterone in circulation is not bound to any protein and is biologically active – this is referred to as ‘free testosterone’ (FT). Bioavailable testosterone refers to the sum of FT and the testosterone bound to serum albumin, since it is bound with low affinity and readily able to dissociate to become available for its biological function. In males elevated levels of testosterone are associated with several conditions such as, early (precocious) puberty, congenital adrenal hyperplasia (CAH), androgen insensitivity syndrome (AIS), steroid use and testicular or adrenal tumours. Whereas the major causes of suppressed levels include Klinefelter’s syndrome, testicular damage, pituitary disorders etc. In females of all ages, elevated testosterone levels can be associated with a variety of virilising conditions including adrenal tumours and polycystic ovarian syndrome (PCOS). These clinical conditions are associated with either a lack or excess of testosterone in circulation (hypoandrogenism or hyperandrogenism). Diagnosis of these disorders involve the quantification of total testosterone (TT) in association with other clinical evidence and laboratory data. However, clinical manifestations of androgen disorders are often associated with normal levels of TT. In such cases, additional information may be gained by the assessment of the biologically active, FT level. Several androgen disorders can be caused by alteration of SHBG production which affects the levels of FT available in serum. Measurement of FT can be considered useful in the diagnosis of several conditions including androgen deficiency in men and androgen excess in women<sup>1</sup>. Assessment of free testosterone levels may prove beneficial<sup>2</sup> and may avoid an incorrect diagnosis of hypogonadism in cases when low concentrations of total testosterone are determined and alterations of SHBG are suspected. There is an observed and well documented circadian variation of testosterone levels in men with circulating concentrations being higher in the morning and declining throughout the day<sup>3</sup>. Testosterone levels also decline in ageing males (andropause) and is associated with loss of muscle and bone mass, leading to osteoporosis, loss of libido, erectile dysfunction. Depression and impaired cognitive function<sup>4</sup>.

## 3. PRINCIPLE OF THE METHOD

The Free Testosterone ELISA is a competitive enzyme immunometric assay (ELISA) where free testosterone (antigen) in the sample competes with the antigenic testosterone conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti- testosterone coated on the microplate (solid phase). After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in the bound fraction reacts with the Substrate (H<sub>2</sub>O<sub>2</sub>) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H<sub>2</sub>SO<sub>4</sub>) is added. The colour intensity is inversely proportional to the free testosterone concentration of in the sample. Free testosterone concentration in the sample is calculated through a calibration curve.

## 4. REAGENTS, MATERIALS AND INSTRUMENTATION

### 4.1. Reagents and materials supplied in the kit

- CAL 0 – 5** Calibrators (6 vial, 1 mL each); Calibrator 0: refer to section 5.
- CONTROL 1 & 2** Controls (2 vials, 1 mL each), Control 1 and 2; refer to section 5. Controls concentration is indicated on the Certificate of Analysis
- ENZ CONJ** Conjugate (1 vial, 15 mL). Testosterone conjugated with Horseradish peroxidase (HRP); ProClin >0.0015% and BSA
- SORB MT** Coated Microplate (1 breakable microplate). Anti-Testosterone antibody adsorbed on microplate
- SUB TMB** TMB Substrate (1 vial, 15 mL). H<sub>2</sub>O<sub>2</sub>-TMB (0.26 g/L) (*avoid any skin contact*). ProClin <0.0015%
- STOP SOLN** Stop Solution (1 vial, 15 mL). Sulphuric acid 0.15M (*avoid any skin contact*)
- WASH SOLN 10x** 10x Conc. Wash Solution (1 vial, 50 mL) 0.2M Phosphate buffer pH 7.4; ProClin >0.0015%

### 4.2. Materials required but not provided

Distilled water

### 4.3. Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

## 5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material. Calibrators and Controls contain Human Serum.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious. Calibrators 1 to 5 and Controls contain BSA (Bovine Serum Albumin)
- Some reagents (calibrator 0, conjugate and wash solution) contain small amounts of ProClin™ 300 (>0,0015%, <0,06%) as preservative. Some reagents (Calibrator 1-5 and TMB Substrate) contain ProClin™ 300 (<0,0015%). Avoid contact with skin or mucosa.
- Classification according to Regulation (EC) No. 1272/2008 [CLP]**

Skin sensitivity, Category 1

Contains: ProClin 300

Warning

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- Some reagents (calibrators and controls) contain small amounts of Sodium Azide ( $\text{NaN}_3$ ) <0.1% which may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through with large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/ $\text{H}_2\text{O}_2$  to directed sunlight, metals or oxidants. Do not freeze the solution.

## 6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

## 7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2-8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

## 8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C	24 hours
Freeze/thaw cycles	1 cycle

## 9. PROCEDURE

### 9.1. Preparation of Calibrators and Controls

Before use, mix for 5 minutes with a rotating mixer.

The Calibrators are ready for use and have the following concentration of Testosterone:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
pg/mL	0	0.2	1.0	4.0	20.0	100.0

The Controls are ready to use.

### 9.2. Preparation of the Conjugate

The conjugate is ready to use. Mix gently, for 5 minutes, on a roller mixer.

### 9.3. Preparation of the Wash Solution

Dilute the content of the vial "Conc. Wash Solution 10X" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

### 9.4. Preparation of Samples

The determination of free testosterone can be performed in human serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

### 9.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C <sub>0</sub> -C <sub>5</sub>	20 µL		
Sample/ Control		20 µL	
Conjugate	100 µL	100 µL	
Incubate for 1 hour at 37°C (± 0.5°C).			
Remove the content from each well; wash the wells 3 times with 300 µL of diluted wash solution			
<b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.			
<b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22÷28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently.			
Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods. The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate. The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate. Demeditec recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range. When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories<sup>5</sup>.

## 11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is required**.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

### Conversion of units

To convert results to SI units:

$$\text{pmol/L} = \text{pg/mL} \times 3.47$$

To convert results to mass units:

$$\text{pg/mL} = \text{pmol/L} \times 0.289$$

## 12. MEASURING RANGE

The assay measuring range (AMR) is 0.34 – 70.0 pg/mL (1.18 – 242.9 pmol/L).

Any value that reads below 0.34 pg/mL (1.18 pmol/L) should be reported as “< 0.34 pg/mL (< 1.18 pmol/L)”. Any value that reads above 70.0 pg/mL (242.9 pmol/L) should be reported as “> 70.0 pg/mL (242.9 pmol/L)”.

## 13. METROLOGY AND TRACEABILITY

The Free Testosterone ELISA has been standardised against internal reference standards (serum matrix) which have been value assigned to another commercially available test method.

#### 14. EXPECTED VALUES

The following ranges were determined using the Free Testosterone ELISA and are provided for information only.

	No. of subjects	Median pg/mL	Reference Interval (pg/mL)
<b>Males</b>			
<b>21 – 49 years</b>	120	14.13	5.01 – 27.78
<b>&gt; 50 years</b>	120	12.75	4.11 – 21.85
<b>Females</b>			
<b>pre-menopausal</b>	120	0.55	< LOQ – 1.70
<b>post-menopausal</b>	120	0.75	< LOQ – 2.34

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

#### 15. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

##### 15.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, “Protocols for Determination of Limits of Detection and Limits of Quantitation” using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.10 pg/mL
Limit of Detection (LoD)	0.20 pg/mL
Limit of Quantitation (LoQ)	0.34 pg/mL

##### 15.2. Trueness

Trueness has been demonstrated through method comparison of the Free Testosterone ELISA to a commercially available assay using native donor samples – refer to section 15.5

##### 15.3. Precision

Precision of the Free Testosterone ELISA was determined by performing a complex precision study.

**Repeatability:** A total of 6 serum samples were assayed in 5 replicates, once per day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (pg/mL)	Within run (Repeatability)	
			SD	CV%
1	75	0.68	0.08	11.3%
2	75	1.62	0.10	5.9%
3	75	5.25	0.19	3.7%
4	75	10.43	0.48	4.6%
5	75	34.99	1.17	3.3%
6	75	68.65	4.62	6.7%

**Reproducibility:** A total of 6 serum samples were assayed in 5 replicates, once per day for 5 days by 3 operators. Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (pg/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	0.75	0.14	18.0%
2	150	1.74	0.24	14.1%
3	150	5.48	0.52	9.4%
4	150	10.72	0.94	8.7%
5	150	37.41	4.83	12.9%
6	150	73.81	10.83	14.7%

**15.4. Linearity**

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For free testosterone concentration by Free Testosterone ELISA, the measurement procedure shows linearity for the interval from 0.27 to 89.83 pg/mL within the allowable deviation of linearity (ADL) of  $\pm 15\%$ .

**15.5. Method comparison**

The Testosterone free ELISA was compared against a commercially available quantitative manual ELISA test, following CLSI EP-9A, "Method Comparison and Bias Estimation Using Patient Samples". A total of 95 samples, selected to represent a wide range of free testosterone concentrations, was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
95	0.94 [0.88 to 0.99]	-0.02 [-0.14 to 0.06]	0.93

**15.6. Analytical Specificity**

The specificity was assessed with the following cross-reactants.

Cross-reactant	Concentration tested (unit)	Mean %Cross reactivity
11-keto-testosterone	100 ng/mL	0.1%
11- $\beta$ -hydroxy-testosterone	10 ng/mL	0.2%
17 $\alpha$ OH-Progesterone	500 ng/mL	0.0%
Aldosterone	3000 ng/mL	0.0%
Androstenedione	100 ng/mL	0.0%
Cortisol	1000 ng/mL	0.0%
Cortisone	1000 ng/mL	0.0%
Danazol	1000 ng/mL	0.0%
Dexamethasone	2000 ng/mL	0.0%
DHEA	1000 ng/mL	0.0%
DHEA-S	10000 ng/mL	0.0%
5 $\alpha$ -dihydrotestosterone (DHT)	500 ng/mL	0.0%
Estradiol	1000 ng/mL	0.0%
Estriol	100 ng/mL	0.0%
Estrone	1000 ng/mL	0.0%
Ethisterone	100 ng/mL	0.0%
Norgestrel	100 ng/mL	0.0%
Prednisone	1000 ng/mL	0.0%
Pregnenolone	5000 ng/mL	0.0%
Progesterone	1000 ng/mL	0.0%
Testosterone propionate	1000 ng/mL	0.0%

The following substances do not interfere with a bias of  $> \pm 15\%$  in the Free Testosterone ELISA when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	200 mg/dL
Total Protein	7 g/dL
Triglyceride	500 mg/dL

**15.7. Serum-plasma study**

The Free Testosterone ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 21 samples (18 native, 3 spiked) to cover the range were evaluated. Linear regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
SST	0.97 [0.92 to 1.02]	-0.02 [-1.01 to 0.98]	0.99
Lithium Heparin	0.99 [0.97 to 1.02]	-0.01 [-0.40 to 1.02]	1.00
Sodium Heparin	0.93 [0.87 to 0.99]	-0.01 [-1.12 to 1.10]	0.99
EDTA	1.03 [0.97 to 1.08]	0.27 [0.74 to 1.28]	0.99

**16. LIMITATIONS OF USE**

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays<sup>6</sup>. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

**17. WASTE MANAGEMENT**

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

**18. BIBLIOGRAPHY**

1. Shea JL, Wong PY, Chen Y. Free testosterone: clinical utility and important analytical aspects of measurement. *Adv Clin Chem.* 2014;63:59-84.
2. Diver MJ. Analytical and physiological factors affecting the interpretation of serum testosterone concentration in men. *Ann Clin Biochem.* 2006 Jan;43(Pt 1):3-12.
3. Brambilla DJ, Matsumoto AM, Araujo AB and McKinlay JB. The Effect of Diurnal Variation on Clinical Measurement of Serum Testosterone and Other Sex Hormone Levels in Men. *J Clin Endocrinol Metab.* 2009 Mar; 94(3): 907–913.
4. Rajfer J. Decreased Testosterone in the Aging Male. *Rev Urol.* 2003;5(suppl 1):S1–S2.
5. Basic QC Practices On-line Course; <http://www.Westgard.com>.
6. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

**19. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT**

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority. The manufacturer can be contacted through their customer service or technical support team. The contact details can be found on the company website:

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7. Shea JL, Wong PY, Chen Y. Free testosterone: clinical utility and important analytical aspects of measurement. *Adv Clin Chem.* 2014;63:59-84.
8. Diver MJ. Analytical and physiological factors affecting the interpretation of serum testosterone concentration in men. *Ann Clin Biochem.* 2006 Jan;43(Pt 1):3-12.
9. Brambilla DJ, Matsumoto AM, Araujo AB and McKinlay JB. The Effect of Diurnal Variation on Clinical Measurement of Serum Testosterone and Other Sex Hormone Levels in Men. *J Clin Endocrinol Metab.* 2009 Mar; 94(3): 907–913.
10. Rajfer J. Decreased Testosterone in the Aging Male. *Rev Urol.* 2003;5(suppl 1):S1–S2.
11. Basic QC Practices On-line Course; <http://www.Westgard.com>.
12. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

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


## SYMBOLS USED WITH DEMEDITEC ASSAYS

Symbol	English	Deutsch	Française	Espanol	Italiano
	European Conformity	CE-Konformitätskennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	utilisation Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di catalogo
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Note warnings and precautions	Warnhinweise und Vorsichtsmaßnahmen beachten	Avertissements et mesures de précaution font attention	Tiene en cuenta advertencias y precauciones	Annoti avvisi e le precauzioni
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
<i>Distributed by</i>	Distributed by	Vertrieb durch	Distribution par	Distribución por	Distribuzione da parte di
	Version	Version	Version	Versión	Versione
	Single-use	Einmalverwendung	À usage unique	Uso único	Uso una volta